

SPECTROPHOTOMETRY AND SPECTROMETRY - CONCEPT AND APPLICATIONS

Renjini A¹, Dani Dileep²

¹ Assistant Professor, Department of ECE, Rajadhani Institute of Engineering and Technology, Kerala, India

² PG scholar, Department of Optoelectronics, University of Kerala, Kerala, India

ABSTRACT

Spectrophotometry is a technique which measures how much a substance absorbs light by measuring the intensity of light, as a beam of light passes through it. The basic principle is that every compound absorbs or transmits light over a certain range of frequencies (wavelengths). If we take a measure of this, it could be used to measure the amount of a known substance in an unknown medium. Spectrophotometry is one of the most useful methods of quantitative analysis in various fields such as chemistry, physics, biochemistry, forensic science, agriculture, material and chemical engineering and clinical applications. A spectrometer gives the detail of exactly which all wavelengths of light are absorbed and which all wavelengths of light are reflected whereas a spectrophotometer gives a comparative measure by measuring the relative intensity of the light absorbed or reflected at a particular wavelength of light. Spectrometric and spectrophotometric techniques can be used for a wide range of elemental and compositional analyses that provide a very high and accurate detection limits in the sub-ppm level of chemical compound which could be made useful in various applications. Some of the spectrophotometers can be used to monitor the separations obtained while doing liquid chromatography and some others could be used for colour selection and printing. Optimum balance of NPK in soil and protein concentration in various substances can also be monitored by the use of spectrometers. Thus there are many key areas where spectrometry and spectrophotometry are used like pharmacology, industries, agriculture, printing, painting and much more of which some of them are detailed in this paper.

Keyword: - *spectrophotometry, spectrometer, chromatography*

1. INTRODUCTION

We know that every compound absorbs, transmits, or reflects light (electromagnetic radiation) over a certain range of wavelength. Spectrophotometry is a measurement of how much a substance absorbs or transmits light [1]. Spectrophotometry is widely used for quantitative analysis in various areas of science and engineering. Any application that deals with substances or materials can use this technique for quantitative as well as qualitative analysis. In biochemistry, for example, it is used to determine enzyme-catalyzed reactions. In clinical applications, it is used to examine blood or tissues for clinical diagnosis [2]. There are also several variants in spectrophotometry such as atomic absorption spectrophotometry and atomic emission spectrophotometry which are used for routine measurements in hospitals, petrochemical industry, food industry, water quality (both purity monitoring and purification) control labs, chemical and biological plants etc. Depending on the range of wavelength of light source, a spectrophotometer can be classified basically into two different types:

- **UV-Visible spectrophotometer:** uses light over the ultraviolet range (185 - 400 nm) and visible range (400-700 nm) of electromagnetic radiation spectrum.
- **IR spectrophotometer:** uses light over the infrared range (700-1300 nm) of electromagnetic radiation spectrum.

Fig 1 shows the UV, IR and visible range of frequencies in the electromagnetic spectrum.




Fig 1: Spectrum showing UV, Visible and IR range of frequencies

In visible spectrophotometry, the absorption or the transmission of a certain substance can be determined by the observed color. For instance, a solution sample that absorbs light over all visible ranges (i.e., transmits none of visible wavelengths) appears to be black in theory. On the other hand, if all visible wavelengths are transmitted (i.e., absorbs nothing), the solution sample appears to be white. If a solution sample absorbs red light (~700 nm), it appears green because green is the complementary color of red [3]. Visible spectrophotometers, in practice, use a prism to attenuate a certain range of wavelength (to filter out other wavelengths) so that the particular beam of light is passed through a sample substance under test [4].

2. STRUCTURE AND MECHANISM

Fig 2 illustrates the basic structure of spectrophotometer. It consists of a source of light, a collimator, a monochromator, a wavelength selector, a cuvette for sample solution, a photo detector, and a digital display or a meter.




Fig 2: Basic structure of spectrophotometer (illustrated by Heesung Shim)

A spectrophotometer, in general, consists of two devices; a spectrometer and a photometer. A spectrometer is a device that produces, typically disperses and measures light. A photometer will usually be a photo detector that measures the intensity of light [4].

- **Spectrometer:** It produces a desired range of wavelength of light. First a collimator (lens) transmits a straight beam of light (photons) that passes through a monochromator (prism) to split it into several component wavelengths (spectrum). Then a wavelength selector (slit) transmits only the desired wavelengths, as shown in Fig 2.
- **Photometer:** After the desired range of wavelength of light passes through the solution of a sample in cuvette, the photometer detects the amount of photons that is absorbed and then sends a signal to a galvanometer or a digital display.

A spectrometer produces a variety of wavelengths because different compounds absorb best at different wavelengths. Once the intensity of light is known after it is made to pass through the cuvette, it could be related to the transmittance (T).

Transmittance is the fraction of light that passes through the sample as shown in Fig 3. This can be calculated using the equation (1) as,

$$\text{Transmittance, } T = I_t/I_0 \text{ ----- (1)}$$

I_t is the light intensity after the beam of light passes through the cuvette and I_0 is the light intensity before the beam of light passes through the cuvette.




Fig 3: Transmittance

As far as IR spectrophotometry is considered, IR radiation does not have enough energy to induce electronic transitions as seen with UV. Absorption of IR is restricted to compounds with small energy differences in the possible vibrational and rotational states. For a molecule to absorb IR, the vibrations or rotations within a molecule must cause a net change in the dipole moment of the molecule. The alternating electrical field of the radiation interacts with fluctuations in the dipole moment of the molecule. If the frequency of the radiation matches the vibrational frequency of the molecule then radiation will be absorbed, causing a change in the amplitude of molecular vibration [5]. These vibrational levels are quantized, and absorption of IR by gases yields a line spectra. However, in liquids or solids, these lines broaden into a continuum due to molecular collisions and other interactions. Transmittance is related to absorption by the expression in equation (2).

$$\text{Absorbance, } A = -\log (T) = -\log (I_t/I_0) \text{ ----- (2)}$$

3. APPLICATIONS

3.1 In Pharmacological Laboratories

In laboratories, UV and Visible spectrophotometers can be used to monitor the separations obtained while doing liquid chromatography. If a mixture is separated in columns, spectrophotometers can be used to detect the different compounds. Using a spectrophotometer is relatively cheap and an easy detector method compared to mass spectrometry detectors. Many illegal drugs contain aromatic groups that produce characteristic UV/Visible spectra. As the absorption peak depends on the nature of the aromatic group(s) in the drug compound, additional spectroscopic information regarding drug derivatives and their metabolites will be obtained.

3.2 In Digital printing

Spectrophotometers are widely used in digital printing technique for colour selection and confirmation. A spectral reflectance curve is obtained by taking the readings at every 10 nm along the visible wavelengths (400 - 700nm). These curves confirm a particular colour to a given specification. The colour management process requires a reference colour data set and softwares for analysing the measured data in addition to a spectrophotometer.

3.3 In Agriculture

In the field of agriculture, production of the crop, water and nutrition management etc. are directly knotted with maintaining the optimum balance of nitrogen, phosphorus, and potassium in the soil for the specific crop concerned with. The correct nutrient levels together with correct soil pH levels also support nutrient availability, organic matter decomposition, and overall soil fertility balance [6]. Spectrophotometer acts as an important tool for monitoring the

nitrogen and phosphorus levels in fertilizers and soil compounds. The molybdenum blue phosphate method in conjugation with UV Visible spectrophotometer is used to determine the amount of phosphate. The added molybdate forms phosphomolybdate complex followed by the reduction of the complex with hydrazine hydrate in aqueous sulphuric acid medium. The amount of phosphate present in the fertilizer is found to be proportional to the color intensity of the reduced phosphomolybdate solution [7]. The presence of nitrogen in the fertilizer can be analyzed rapidly by adding several drops of concentrated FeSO_4 solution. The nitrogen levels are quantified by a spectrophotometer based on the hue of the solution and the color density measurements.

3.4 Estimation of protein concentration

Due to the presence of tyrosine and tryptophan, the proteins absorb UV light at a maximum wavelength of 280 nm. The tyrosine and tryptophan content varies with different proteins and hence the absorption peak. The absorption peak at 280 nm is hence used as a rapid and a very sensitive measure of protein concentration [8]. With the amount of known absorbance the unknown concentration of the sample can be found by using Beer-Lambert Law (also known as Beer's Law) states that there is a linear relationship between the absorbance and the concentration of a sample. Beer's Law could be expressed in equation form as given by equation (3).

$$A = \epsilon lc \text{ ----- (3)}$$

where A is the measure of absorbance, ϵ is the molar extinction coefficient or molar absorptivity (or absorption coefficient), l is the path length, and c is the concentration [4]. Even though we could obtain accurate readings from a spectrometer/ spectrophotometer, there are some common limitations. Dust or grim coated over the mirrors in a double-beam UV-Visible spectrometer may degrade the device's performance. Also the stray light caused by faulty equipment design could decrease an instrument's linearity range and reduce the absorbency of the substance it measures. Another factor is that the electronic components in the spectrometer or sample source may generate noise that decreases measurement accuracy and reduce the device sensitivity.

4. CONCLUSION

To conclude, spectroscopic and spectrometric solutions offer very high accuracy and lend a helping hand in eliminating human error with high-speed wavelength scanning of samples across the UV and Visible spectrum. Other than these analyses, techniques that use X-ray Fluorescence technology provide an accurate analytical method finding out the chemical composition of many substances. Electron spectroscopy is yet another analytical technique used to study the electronic structure and its dynamics in atoms and molecules. These analytical methods are non-destructive and reliable, require no or very little, sample preparation and are suitable for solid, liquid and powdered samples. They can be used for a wide range of elemental analysis and provide detection limits at the sub-ppm level measuring concentrations of near hundred percentage easily and simultaneously.

5. REFERENCES

- [1]. https://chem.libretexts.org/Core/Physical_and_Theoretical_Chemistry/Kinetics/Reaction_Rates/Experimental_Determination_of_Kinetics/Spectrophotometry (last accessed on 24th February 2017).
- [2]. <https://www.researchgate.net/post/> (accessed on 15th February 2017).
- [3]. Gore, Michael, Spectrophotometry and Spectrofluorimetry, Newyork: Oxford University Press, 2000.
- [4]. <https://chem.libretexts.org/> (last accessed on 19th February 2017).
- [5]. N Raaman, "Phytochemical Techniques", 2006, ISBN 81-89422-30-1.
- [6]. <https://in.hach.com/> (accessed on 18th February 2017).
- [7]. Samjhana Pradhan, Megh Raj pokhrel, "Spectrophotometric determination of phosphate in in sugarcane juice, fertilizer ,detergent and water samples by molybdenum blue method", Scientific world, vol.11, July 2013.
- [8]. Whitaker, Granum, "An absolute method for protein determination based on difference in absorbance at 235 and 280 nm", *Analyt. Biochem.*, 1980.

BIOGRAPHIES

	<p>Ms Renjini A is currently working as Assistant Professor in Department of ECE, Rajadhani Institute of Engineering and Technology, Trivandrum. She is having a close association with this esteemed institution for seven years. Her areas of interest include optical and microwave communication, biomedical engineering as well as instrumentation and nanotechnology. She has several national and international publications including an IEEE journal to her credit. Her paper titled “Three Dimensional Holographic Data Storage” bagged first prize in the Sixth International Conference on Computing and Communication Technology held at IAJJCE, Marthandam on March 9th, 2017.</p>
	<p>Ms Dani Dileep is a PG scholar at Department of Optoelectronics, University of Kerala, Karyavattom. Recently she had presented a paper entitled “Study of Carbon Nanoparticles from Diesel Engine” in the National conference on Recent Trends in Engineering and Technology held at Trinity College of Engineering, Trivandrum on February 20,2017.</p>